


LVA30

CATALOGUE 2024-25



HSN Code : 90278910



## Low Viscosity Adaptor


LVA30

- Spindle range - 1~2K for LMDV-60
- Spindle range - 0.6~6K for LMDV100
- Spindle range - 0.3~6K for LMDV200
- K=1000



### SALIENT FEATURES

- Ideal for Low Viscosity Materials
- 20 to 30 ml sample is needed for test in Low Viscosity Adaptor (LVA30)
- Spindle Code.L0 is used with Low Viscosity Adaptor
- Stainless steel parts for easy to clean



**Note :** LVA30 is used for 1 to 50m.Pas viscosity measuring range.

**STANDARD ACCESSORIES :** Main Instrument, Joint rod, Hook rod, Frame, Spindle Code. L0, L0 Pipe