

SSA15

CATALOGUE 2024-25



## Small Sample Adaptor

SSA15

The Small Sample Adaptor, consisting of a coupling nut, an extension link, a spindle set (spindle#21,#27,#28 and #29), a sample chamber and an adjustable base

### SALIENT FEATURES


- Small Sample Adaptor is recommended for small sample volume.
- Small Sample Adaptor widely used for Digital Viscometer
- 10 to 20 ml sample is needed for test in Small Sample Adaptor.
- Sample temperature can be controlled by using Refrigerated Circulating Bath (Optional)
- The Sample Chamber fit into a water jacket.

### Spindle Range Table


HSN Code : 90278910

Spindle Model	21	27	28	29
LMDV60	17~ 16.7K	83~ 83.3K	167~ 166.7K	333~ 333.3K
LMDV100	10~50K	50~250K	100~500K	200~1M
LMDV200	5~50K	25~250K	50~500K	100~1M

K=1000 M=1000000



### Spindle set



\* Small Sample Adaptor is an Optional Accessory with Digital Viscometer

Note : Small Sample Adaptor order along with Viscometer only. It **cannot** be attached with the main instrument later on.

**STANDARD ACCESSORIES :** Main Instrument, Adjustable Base, Sample Chamber, Extension & Coupling Nut, Spindle Set

★ Specifications are subject to change without notice.